Gebrauchsanleitung TPP – Clipmax / Instructions for Use TPP – Clipmax

Anleitung

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Instructions

Clipmax is a tissue culture flask 10 cm² with the standard microscope slide Clipmax (25 x 75 mm). It is intended for cell culture, staining, and immunofluorescence as it does not interfere with fluorescent labels-directly on the slide. It is useful for viral and mycoplasma testing, chromosome studies toxicity tests and immunocytology.

The growth area is opto-mechanically treated for an optimal cell adhesion and proliferation.

After the cell culture the medium chamber with the biocompatible seal ring can easily be detached from the slide. Cells may then be fixed, stained; cover slipped, and filed with the usual convenience and can be used directly for microscopic view. Clipmax features an excellent chemical resistance and thus allows cell fixation and the staining of organic solvents for EM and other techniques.

The system allows a constant sterile gas exchange through the filter screw cap.

- Chamber separation by clip open, no separation tool required
- Sealant: Biocompatible acrylic adhesive remains with media chamber after detachment
- Open product packaging
- Inspect for intact bag and product.
- Open filter screw-cap
- Fill medium: Introduce medium or inoculate cells with a pipette alongside the wall, or direct injection on the growth area of the Clipmax.
- Close flask with filter screw cap
- Incubate at 37 °C in an incubator, little nubs on the underside of the slide protect from scratches during mechanical movements
- Aspirate medium with pipette
- Detach medium chamber
- Fix, stain and harvest cells following protocol

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Advice and Recommendation

For homogenous cultivation of the cells it is most important to level the incubator and the continuous supply of medium. Sources of vibrations must be eliminated to avoid circular growth pattern. The necessary gas exchange takes place through the screw cap.

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Materialien

Medien-Kammer PS
 Clipmax COP
 Kappe PE
 Membrane PTFE,
 Porengrösse 0.22 µm

Technical Data

Materials

Medium chamber PS
 Slide COP
 Screw cap PE
 Membrane PTFE,
 pore size 0.22 µm

Abmessungen	Measurement	70010
Kappe	Cap	Filter
Membrane µm	Membrane µm	0.22
Höhe mm	Heights mm	32
Breite mm	Width mm	33
Länge mm	93	92
Wachstums-fläche cm ²	Growth area cm ²	10
Empf Volumen ml	Rec. volume ml	3.5
Max Volumen ml	Max. volume ml	5
Clipmax Breite x Länge mm	Clipmax Width x Length mm	25 x 75